

Anti-Lipoprotein lipase LPL Rabbit Monoclonal Antibody

Catalog Number: M00316


Overview

Product Name	Anti-Lipoprotein lipase LPL Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Lipoprotein lipase LPL Rabbit Monoclonal Antibody catalog # M00316. Tested in WB, IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal ABDE-12
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P06858


Technical Details

Immunogen	A synthesized peptide derived from human Lipoprotein lipase
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:1000-1:5000 IHC 1:50-1:200


Anti-Lipoprotein lipase LPL Rabbit Monoclonal Antibody (M00316) Images




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Lipoprotein lipase expression in Human fetal liver lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

2 Publications Citing This Product

1. PubMed ID: 30456017, Supplemental Information

2. PubMed ID: 21762526, Liu Y, Wang Zb, Yin Wd, Li Qk, Cai Mb, Yu J, Li Hg, Zhang C, Zu Xh. Lipids Health Dis. 2011 Jul 16;10:117. Doi: 10.1186/1476-511X-10-117. Preventive Effect Of Ibroloipim On Suppressing Lipid Accumulation And Increasing Lipoprotein Lipase In The Kid...

Visit bosterbio.com/anti-lipoprotein-lipase-rabbit-monoclonal-antibody-m00316-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Lipoprotein lipase LPL Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.