

Anti-Thyroid Peroxidase Rabbit Monoclonal Antibody

Catalog Number: M00320-2


Overview

Product Name	Anti-Thyroid Peroxidase Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Thyroid Peroxidase Rabbit Monoclonal Antibody catalog # M00320-2. Tested in WB, IHC, IP applications. This antibody reacts with Human.
Application	IP, IHC, WB
Clonality	Monoclonal 24T05
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07202

Technical Details

Immunogen	A synthesized peptide derived from human Thyroid Peroxidase
Isotype	IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB 1:1000-1:5000 IHC 1:100-1:500 IP 1:50</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-Thyroid Peroxidase Rabbit Monoclonal Antibody (M00320-2) Images



Western blot analysis of Thyroid Peroxidase expression in Human thyroid lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Thyroid Peroxidase Rabbit Monoclonal Antibody