

## Anti-Phospho-JAK1 (Y1034 + Y1035) Rabbit Monoclonal Antibody

Catalog Number: M00330Y1034Y1035


### Overview

Product Name	Anti-Phospho-JAK1 (Y1034 + Y1035) Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Phospho-JAK1 (Y1034 + Y1035) Rabbit Monoclonal Antibody catalog # M00330Y1034Y1035. Tested in WB, IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal 31J16
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P23458

### Technical Details

Immunogen	A synthesized peptide derived from human Phospho-JAK1 (Y1034 + Y1035)
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-100

## Anti-Phospho-JAK1 (Y1034 + Y1035) Rabbit Monoclonal Antibody (M00330Y1034Y1035) Images



Western blot analysis of Phospho-JAK1 (Y1034 + Y1035) expression in Ramos treated with pervanadate cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-JAK1 (Y1034 + Y1035) Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.