

Anti-MTHFR Monoclonal Antibody

Catalog Number: M00336


Overview

Product Name	Anti-MTHFR Monoclonal Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-MTHFR Monoclonal Antibody catalog # M00336. Tested in ELISA, IHC, WB applications. This antibody reacts with Human, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 5D3
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P42898


Technical Details

Immunogen	Purified recombinant fragment of human MTHFR expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-MTHFR Monoclonal Antibody (M00336) Images



Immunohistochemistry analysis of paraffin-embedded lung cancer tissues with DAB staining using MTHFR Monoclonal Antibody.



Western Blot analysis using MTHFR Monoclonal Antibody against HEK293 (1) and MTHFR-hIgGfc transfected HEK293 (2) cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MTHFR Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.