

Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (tisotumab)

Catalog Number: M00342-3


Overview

Product Name	Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (tisotumab)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (tisotumab) (Catalog # M00342-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.332EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P13726


Technical Details

Isotype	IgG1
Form	Liquid

Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (tisotumab) (M00342-3) Images



Anti-F3/Factor III/Tissue Factor/CD142 Reference Antibody (tisotumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-F3/Factor III/Tissue Factor/CD142 Reference Antibody (tisotumab) is more than 99.08%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (tisotumab)

For Research Use Only. Not for use in diagnostic procedures.