

Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (Chugai patent anti-TF)

Catalog Number: M00342-5


Overview

Product Name	Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (Chugai patent anti-TF)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (Chugai patent anti-TF) (Catalog # M00342-5). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.334EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P13726


Technical Details

Isotype	IgG4SP
Form	Liquid

Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (Chugai patent anti-TF) (M00342-5) Images



Anti-F3/Factor III/Tissue Factor/CD142 Reference Antibody (Chugai patent anti-TF) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-F3/Factor III/Tissue Factor/CD142 Reference Antibody (Chugai patent anti-TF) is more than 95%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-F3 / Factor III / Tissue Factor / CD142 Reference Antibody (Chugai patent anti-TF)

For Research Use Only. Not for use in diagnostic procedures.