

Anti-CEA (CD66e) CEACAM5 Rabbit Monoclonal Antibody

Catalog Number: M00356-1


Overview

Product Name	Anti-CEA (CD66e) CEACAM5 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-CEA (CD66e) CEACAM5 Rabbit Monoclonal Antibody catalog # M00356-1. Tested in WB, IHC, ICC/IF applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Monoclonal AAOA-3
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P06731


Technical Details

Immunogen	A synthesized peptide derived from human CEA (CD66e)
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IHC 1:50-1:200 ICC/IF 1:50-1:200


Anti-CEA (CD66e) CEACAM5 Rabbit Monoclonal Antibody (M00356-1) Images




Western blot analysis of CEA(CD66e) expression in Human colon cancer lysate.



Immunofluorescent analysis using the Antibody at 1:50 dilution.




Immunohistochemical analysis of paraffin-embedded human colon carcinoma, using CEA(CD66e) Antibody.



Immunofluorescent analysis using the Antibody at 1:150 dilution.

Immunofluorescent analysis of MCF-7 cells, using CEA(CD66e) Antibody.



Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CEA (CD66e) CEACAM5 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.