

Anti-CEACAM5 / CEA / CD66e Reference Antibody (Immunomedics patent anti-CEACAM5 (Class III))

Catalog Number: M00356-8


Overview

Product Name	Anti-CEACAM5 / CEA / CD66e Reference Antibody (Immunomedics patent anti-CEACAM5 (Class III))
Reactive Species	Human
Description	Boster Bio Anti-CEACAM5 / CEA / CD66e Reference Antibody (Immunomedics patent anti-CEACAM5 (Class III)) (Catalog # M00356-8). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.346EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P06731


Technical Details

Isotype	IgG1
Form	Liquid

Anti-CEACAM5 / CEA / CD66e Reference Antibody (Immunomedics patent anti-CEACAM5 (Class III)) (M00356-8) Images



Anti-CEACAM5/CEA/CD66e Reference Antibody (Immunomedics patent anti-CEACAM5 (Class III)) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-CEACAM5/CEA/CD66e Reference Antibody (Immunomedics patent anti-CEACAM5 (Class III)) is more than 95%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CEACAM5 / CEA / CD66e Reference Antibody (Immunomedics patent anti-CEACAM5 (Class III))

For Research Use Only. Not for use in diagnostic procedures.