

## Anti-Thyroglobulin TG Monoclonal Antibody

Catalog Number: M00359-4

### About TG

C3 plays a central role in the activation of the complement system. Its processing by C3 convertase is the central reaction in both classical and alternative complement pathways. After activation C3b can bind covalently, via its reactive thioester, to cell surface carbohydrates or immune aggregates.

### Overview

Product Name	Anti-Thyroglobulin TG Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Thyroglobulin TG Monoclonal Antibody catalog # M00359-4. Tested in WB, IHC, ICC/IF, IP applications. This antibody reacts with Human.
Application	IP, IF, IHC, ICC, WB
Clonality	Monoclonal AFEI-20
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01266

### Technical Details

Immunogen	A synthesized peptide derived from human Thyroglobulin Precursor of the iodinated thyroid hormones thyroxine (T4) and triiodothyronine (T3) .
Isotype	Rabbit IgG
Form	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Affinity-chromatography
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:


WB 1:1000-1:5000

IHC 1:50-1:200

ICC/IF 1:50-1:200

IP 1:50

## Anti-Thyroglobulin TG Monoclonal Antibody (M00359-4) Images



Western blot analysis of Thyroglobulin expression in Human thyroid lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Thyroglobulin TG Monoclonal Antibody