

## Anti-PDGFR alpha Rabbit Monoclonal Antibody

Catalog Number: M00366


### Overview

Product Name	Anti-PDGFR alpha Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PDGFR alpha Rabbit Monoclonal Antibody catalog # M00366. Tested in WB application. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Monoclonal AAOE-16
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P16234

### Technical Details

Immunogen	A synthesized peptide derived from human PDGFR alpha
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

## Anti-PDGFR alpha Rabbit Monoclonal Antibody (M00366) Images



Western blot analysis of PDGFR alpha/PDGFR alpha using anti-PDGFR alpha/PDGFR alpha antibody (M00366). Electrophoresis was performed on a 8% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human Caco-2 whole cell lysates, Lane 2: human SH-SY5Y whole cell lysates, Lane 3: rat ovary tissue lysates, Lane 4: rat brain tissue lysates, Lane 5: mouse ovary tissue lysates, Lane 6: mouse brain tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-PDGFR alpha/PDGFR alpha antigen affinity purified monoclonal antibody (M00366) at 1:500 overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for PDGFR alpha/PDGFR alpha at approximately 180 kDa. The expected band size for PDGFR alpha/PDGFR alpha is at 123 kDa.

## 2 Publications Citing This Product

1. PubMed ID: 10.1016/j.anireprosci.2019.05.006, Effect of maternal feed restriction in dairy goats at different stages of gestation on skeletal muscle development and energy metabolism of kids at the time of births
2. PubMed ID: 29140299, PDGFR $\alpha$  regulated by miR-34a and FoxO1 promotes adipogenesis in porcine intramuscular preadipocytes through Erk signaling pathway

Visit [bosterbio.com/anti-pdgfr-alpha-rabbit-monoclonal-antibody-m00366-boster.html](https://www.bosterbio.com/anti-pdgfr-alpha-rabbit-monoclonal-antibody-m00366-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PDGFR alpha Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.