

Anti-ANGPT2 Reference Antibody (Roche patent anti-ANG-2)

Catalog Number: M00370-5

Overview

Product Name	Anti-ANGPT2 Reference Antibody (Roche patent anti-ANG-2)
Reactive Species	Human
Description	Boster Bio Anti-ANGPT2 Reference Antibody (Roche patent anti-ANG-2) (Catalog # M00370-5). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.362EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	O15123

Technical Details

Isotype	IgG1
Form	Liquid

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ANGPT2 Reference Antibody (Roche patent anti-ANG-2)

For Research Use Only. Not for use in diagnostic procedures.