

Anti-Myeloperoxidase Rabbit Monoclonal Antibody

Catalog Number: M00372-3


Overview

Product Name	Anti-Myeloperoxidase Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of Myeloperoxidase (MPO) protein. Validated for research with IHC, WB, ELISA applications. Anti Myeloperoxidase (MPO) is reactive for Human, Mouse, Rat samples
Application	ELISA, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P05164

Technical Details

Isotype	IgG, Kappa
Purification	Protein A

Anti-Myeloperoxidase Rabbit Monoclonal Antibody (M00372-3) Images



Human spleen tissue was stained with anti-Myeloperoxidase(MPO) rabbit Antibody

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Myeloperoxidase Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.