

Anti-IL-2 Monoclonal Antibody

Catalog Number: M00387


Overview

Product Name	Anti-IL-2 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-IL-2 Monoclonal Antibody catalog # M00387. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 8D5
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P60568

Technical Details

Immunogen	Purified recombinant fragment of IL-2 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-IL-2 Monoclonal Antibody (M00387) Images



Western blotting validation for Anti-IL-2 Monoclonal Antibody M00387 Western Blot (WB) analysis using IL-2 Monoclonal Antibody against full-length IL-2 recombinant protein with Trx tag (1) and full-length IL-2-hlgGfc transfected HEK293 cell lysate(2). Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-2 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.