

Anti-HNF-4-alpha Rabbit Monoclonal Antibody

Catalog Number: M00389


Overview

Product Name	Anti-HNF-4-alpha Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-HNF-4-alpha Rabbit Monoclonal Antibody catalog # M00389. Tested in WB, IHC, ICC/IF, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Monoclonal GBA-8
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P41235


Technical Details

Immunogen	A synthesized peptide derived from human HNF-4-alpha
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IHC 1:50-1:200 ICC/IF 1:50-1:200 FC 1:50

Anti-HNF-4-alpha Rabbit Monoclonal Antibody (M00389) Images



Western blot analysis of HNF-4-alpha expression in A549 cell lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using HNF-4-alpha Antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HNF-4-alpha Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.