

Anti-Retinoic Acid Receptor alpha RARA Rabbit Monoclonal Antibody

Catalog Number: M00392


Overview

Product Name	Anti-Retinoic Acid Receptor alpha RARA Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Retinoic Acid Receptor alpha RARA Rabbit Monoclonal Antibody catalog # M00392. Tested in WB, Flow Cytometry applications. This antibody reacts with Human, Mouse.
Application	Flow Cytometry, WB
Clonality	Monoclonal GBC-18
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P10276


Technical Details

Immunogen	A synthesized peptide derived from human Retinoic Acid Receptor alpha
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 FC 1:20


Anti-Retinoic Acid Receptor alpha RARA Rabbit Monoclonal Antibody (M00392) Images



Western blot analysis of Retinoic Acid Receptor alpha expression in MCF-7 cell lysate.



All lanes use the Antibody at 1:1000 dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1000 dilution for 1 hour at room temperature.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Retinoic Acid Receptor alpha RARA Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.