

Anti-CK-5/CK-6 (Cytokeratin 5/ Cytokeratin 6) Rabbit Monoclonal Antibody, Clone#RM341

Catalog Number: M00398-4


Overview

Product Name	Anti-CK-5/CK-6 (Cytokeratin 5/ Cytokeratin 6) Rabbit Monoclonal Antibody, Clone#RM341
Reactive Species	Human
Description	Boster Bio Anti-CK-5/CK-6 (Cytokeratin 5/ Cytokeratin 6) Rabbit Monoclonal Antibody, Clone#RM341 (Catalog # M00398-4). Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal RM341
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P13647


Technical Details

Immunogen	A peptide corresponding to the C-terminus of Human Cytokeratin 5
Cross Reactivity	This antibody reacts to both human CK-5 (Cytokeratin 5) and CK-6 (Cytokeratin 6).
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	Immunohistochemistry (IHC): 1:500 - 1:1000 dilution WB: 1:2000-1:4000 dilution.

Anti-CK-5/CK-6 (Cytokeratin 5/ Cytokeratin 6) Rabbit Monoclonal Antibody, Clone#RM341 (M00398-4) Images



Western Blotting result Western Blot of A431 cell lysates using RM341 using anti-CK-5/CK-6 rabbit monoclonal antibody (Clone RM341) at a 1:2000 dilution.



IHC result Immunohistochemical staining of formalin fixed and paraffin embedded human tonsil tissue section using anti-CK-5/CK-6 rabbit monoclonal antibody (Clone RM341) at a 1:1000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CK-5/CK-6 (Cytokeratin 5/ Cytokeratin 6) Rabbit Monoclonal Antibody, Clone#RM341

For Research Use Only. Not for use in diagnostic procedures.