

## Anti-TBX21 Rabbit Monoclonal Antibody

Catalog Number: M00404-1

### Overview


Product Name	Anti-TBX21 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-TBX21 Rabbit Monoclonal Antibody catalog # M00404-1. Tested in IHC application. This antibody reacts with Human.
Application	IHC
Clonality	Monoclonal 20T77
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UL17

### Technical Details

Immunogen	A synthesized peptide derived from human T-bet / Tbx21
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	IHC 1:50-200

## Anti-TBX21 Rabbit Monoclonal Antibody (M00404-1) Images

---



Immunohistochemical analysis of paraffin-embedded human spleen, using T-bet/Tbx21 Antibody.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TBX21 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.