

Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (tigatuzumab)

Catalog Number: M00410-6


Overview

Product Name	Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (tigatuzumab)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (tigatuzumab) (Catalog # M00410-6). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.372EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	O14763


Technical Details

Isotype	IgG1
Form	Liquid

Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (tigatuzumab) (M00410-6) Images



Anti-TNFRSF10B/TRAILR2/CD262 Reference Antibody (tigatuzumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-TNFRSF10B/TRAILR2/CD262 Reference Antibody (tigatuzumab) is more than 95%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFRSF10B / TRAILR2 / CD262 Reference Antibody (tigatuzumab)

For Research Use Only. Not for use in diagnostic procedures.