

## Anti-Blimp-1 PRDM1 Monoclonal Antibody

Catalog Number: M00412


### Overview

Product Name	Anti-Blimp-1 PRDM1 Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Blimp-1 PRDM1 Monoclonal Antibody catalog # M00412. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Monoclonal 50000000
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O75626

### Technical Details

Immunogen	Purified recombinant fragment of human Blimp-1 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-Blimp-1 PRDM1 Monoclonal Antibody (M00412) Images



Western Blot analysis using Blimp-1 Monoclonal Antibody against Raji (1, 2), L1210 (3) and TPH-1 (4) cell lysate.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Blimp-1 PRDM1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.