

## Anti-CTLA-8 / IL-17a Antibody (SY18-VHH-11)

Catalog Number: M00421-6


### Overview

Product Name	Anti-CTLA-8 / IL-17a Antibody (SY18-VHH-11)
Reactive Species	Human
Description	Boster Bio Anti-CTLA-8 / IL-17a Antibody (SY18-VHH-11) (Catalog # M00421-6). Tested in Flow Cytometry, ELISA, FTA applications. This antibody reacts with Human. Endotoxin: < 0.001EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, FTA
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q16552


### Technical Details

Isotype	IgG1
Form	Liquid
Concentration	1mg/ml

## Anti-CTLA-8 / IL-17a Antibody (SY18-VHH-11) (M00421-6) Images



Anti-CTLA-8/IL-17a Antibody (SY18-VHH-11) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-CTLA-8/IL-17a Antibody (SY18-VHH-11) is more than 95%

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CTLA-8 / IL-17a Antibody (SY18-VHH-11)

For Research Use Only. Not for use in diagnostic procedures.