

## Anti-AIM2 Antibody

Catalog Number: M00443

### About AIM2

Tumor suppressor which may act by repressing NF-kappa-B transcriptional activity.


### Overview

Product Name	Anti-AIM2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-AIM2 Antibody (Catalog # M00443). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Monoclonal 812CT9.2.2.1
Formulation	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Mouse
Uniprot ID	O14862


### Technical Details

Immunogen	This AIM2 Monoclonal antibody is generated from mice immunized with a KLH conjugated synthetic peptide selected from the 10-40 region of human AIM2.
Predicted Reactive Species	Bovine, Mouse
Isotype	IgM
Purification	This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Suggested Dilutions	WB: 1:1000

## Anti-AIM2 Antibody (M00443) Images



Western blot analysis of lysate from NCI-H292 cell line, using AIM2 Antibody. AIM2 Antibody was diluted at 1:1000. A goat anti-mouse IgG H&L (HRP) at 1:3000 dilution was used as the secondary antibody. Lysate at 20ug.



Western blot analysis of lysate from NCI-H292 cell line, using AIM2 Antibody. AIM2 Antibody was diluted at 1:1000. A goat anti-mouse IgG H&L (HRP) at 1:3000 dilution was used as the secondary antibody. Lysate at 20ug.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-AIM2 Antibody

For Research Use Only. Not for use in diagnostic procedures.