

Anti-NPM1 monoclonal Antibody

Catalog Number: M00450


Overview

Product Name	Anti-NPM1 monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-NPM1 monoclonal Antibody catalog # M00450. Tested in WB,ICC applications. This antibody reacts with Human,Mouse.
Application	ICC, WB
Clonality	Monoclonal 40000000
Formulation	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P06748


Technical Details

Immunogen	Recombinant full length Human NPM1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000 ICC: 1:100-200

Anti-NPM1 monoclonal Antibody (M00450) Images



Western blot detection of NPM1 antibody in HeLa, Jurkat&3T3 cell lysates using NPM1 antibody (1:1000 diluted) .



Immunocytochemistry of HeLa cells using anti- NPM1 antibody diluted 1:200

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NPM1 monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.