

Anti-Cytochrome P450 2C9 Rabbit Monoclonal Antibody

Catalog Number: M00465-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-Cytochrome P450 2C9 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Cytochrome P450 2C9 Rabbit Monoclonal Antibody catalog # M00465-1. Tested in WB, ICC/IF, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, ICC, WB
Clonality	Monoclonal 20C30
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P11712

Technical Details

Immunogen	A synthesized peptide derived from human Cytochrome P450 2C9
Isotype	IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000

	ICC/IF 1:50-200 FC 1:500
--	-----------------------------

Anti-Cytochrome P450 2C9 Rabbit Monoclonal Antibody (M00465-1) Images



Western blot analysis of Cytochrome P450 2C9 expression in HepG2 cell lysate.

Anti-Cytochrome P450 2C9 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.