

## Anti-Integrin $\alpha$ 4/ITGA4/CD49d Reference Antibody (natalizumab)


Catalog Number: M00468-2

### Overview


Product Name	Anti-Integrin $\alpha$ 4/ITGA4/CD49d Reference Antibody (natalizumab)
Reactive Species	Human
Description	Boster Bio Anti-Integrin $\alpha$ 4/ITGA4/CD49d Reference Antibody (natalizumab) (Catalog # M00468-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.012EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P13612

### Technical Details

## Anti-Integrin $\alpha$ 4/ITGA4/CD49d Reference Antibody (natalizumab) (M00468-2) Images



Anti-Integrin  $\alpha$ 4/ITGA4/CD49d Reference Antibody (natalizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-Integrin  $\alpha$ 4/ITGA4/CD49d Reference Antibody (natalizumab) is more than 98.72%

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin  $\alpha$ 4/ITGA4/CD49d Reference Antibody (natalizumab)

For Research Use Only. Not for use in diagnostic procedures.