

Anti-Integrin a4b7 (ITGA4 & ITGB7) Reference Antibody (vedolizumab)


Catalog Number: M00468-3

Overview


Product Name	Anti-Integrin a4b7 (ITGA4 & ITGB7) Reference Antibody (vedolizumab)
Reactive Species	Human
Description	Boster Bio Anti-Integrin a4b7 (ITGA4 & ITGB7) Reference Antibody (vedolizumab) (Catalog # M00468-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.018EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P13612/Q13683

Technical Details

Anti-Integrin a4b7 (ITGA4 & ITGB7) Reference Antibody (vedolizumab) (M00468-3) Images



Anti-Integrin a4b7 (ITGA4 & ITGB7) Reference Antibody (vedolizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-Integrin a4b7 (ITGA4 & ITGB7) Reference Antibody (vedolizumab) is more than 100%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Integrin a4b7 (ITGA4 & ITGB7) Reference Antibody (vedolizumab)

For Research Use Only. Not for use in diagnostic procedures.