

Anti-Glucocorticoid Receptor Rabbit Monoclonal Antibody

Catalog Number: M00503-4

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Glucocorticoid Receptor Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	This Recombinant Rabbit Monoclonal Antibody detects endogenous levels of Glucocorticoid Receptor protein. Validated for research with WB,IHC,IF,IP,ELISA applications. Anti Glucocorticoid Receptor is reactive for Human, Mouse, Rat samples
Application	ELISA, IP, IF, IHC, WB
Clonality	Monoclonal
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% stabilizing protein This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C/1 year
Host	Rabbit
Uniprot ID	P04150


Technical Details

Isotype	IgG, Kappa
Purification	Protein A
Suggested Dilutions	IHC 1:1000-1:4000 WB 1:1000-1:5000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200


Anti-Glucocorticoid Receptor Rabbit Monoclonal Antibody (M00503-4) Images




Rat liver was stained with anti-Glucocorticoid Receptor rabbit antibody



Human liver was stained with anti-Glucocorticoid Receptor rabbit antibody



Mouse liver was stained with anti-Glucocorticoid Receptor rabbit antibody



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Glucocorticoid Receptor antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: U-251 MG Lane 2: A549 Lane 3: RAW264.7
Predicted band size: 85kDa Observed band size: 94kDa

Anti-Glucocorticoid Receptor Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.