

## Anti-Protein delta homolog 1 DLK1 Monoclonal Antibody

Catalog Number: M00513


### Overview

Product Name	Anti-Protein delta homolog 1 DLK1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Protein delta homolog 1 DLK1 Monoclonal Antibody catalog # M00513. Tested in ELISA, Flow Cytometry, IF, WB applications. This antibody reacts with Human.
Application	ELISA, Flow Cytometry, IF, WB
Clonality	Monoclonal 3A10
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P80370


### Technical Details

Immunogen	Purified recombinant fragment of human DLK1 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 FC 1:200-1:400 ELISA 1:10000


## Anti-Protein delta homolog 1 DLK1 Monoclonal Antibody (M00513) Images




Western blotting validation for Anti-Protein delta homolog 1 DLK1 Monoclonal Antibody M00513 Western Blot (WB) analysis using DLK1 Monoclonal Antibody against human DLK1 (AA: 174-349) recombinant protein. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Immunofluorescent staining data of DLK1 using Anti-Protein delta homolog 1 DLK1 Monoclonal Antibody (M00513). Immunofluorescence (IF) analysis of U251 cells using DLK1 Monoclonal Antibody (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow Cytometry validation of DLK1 using Anti-Protein delta homolog 1 DLK1 Monoclonal Antibody (M00513). Flow cytometric (FCM) analysis of NIH/3T3 cells using DLK1 Monoclonal Antibody (green) and negative control (red).



ELISA validation of DLK1 using Anti-Protein delta homolog 1 DLK1 Monoclonal Antibody (M00513). ELISA analysis of DLK1 antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein delta homolog 1 DLK1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.