

Anti-Stat6 Monoclonal Antibody

Catalog Number: M00523-2


Overview

Product Name	Anti-Stat6 Monoclonal Antibody
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-Stat6 Monoclonal Antibody catalog # M00523-2. Tested in ELISA, IF, WB applications. This antibody reacts with Human, Monkey, Mouse, Rat.
Application	ELISA, IF, WB
Clonality	Monoclonal 7D3
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P42226


Technical Details

Immunogen	Purified recombinant fragment of human Stat6 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000

Anti-Stat6 Monoclonal Antibody (M00523-2) Images



Western Blot analysis using Stat6 Monoclonal Antibody against HEK293 (1), NIH/3T3 (2), MCF-7 (3), Raw246.7 (4) and PC-12 (5) cell lysate.



Immunofluorescence analysis of HeLa cells using Stat6 Monoclonal Antibody (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Stat6 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.