

Anti-Ikaros Monoclonal Antibody

Catalog Number: M00531-1


Overview

Product Name	Anti-Ikaros Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Ikaros Monoclonal Antibody catalog # M00531-1. Tested in WB, IHC, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IHC, WB
Clonality	Monoclonal AEFA-9
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13422

Technical Details

Immunogen	A synthesized peptide derived from human Ikaros Transcription regulator of hematopoietic cell differentiation (PubMed:17934067) . Binds gamma-satellite DNA (PubMed:17135265, PubMed:19141594) . Plays a role in the development of lymphocytes, B- and T-cells.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 IHC 1:50-200 FC 1:60

Anti-Ikaros Monoclonal Antibody (M00531-1) Images



Western blot analysis of Ikaros expression in MOLT4 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ikaros Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.