

## Anti-ErbB (HER3) Rabbit Monoclonal Antibody

Catalog Number: M00539


### Overview

Product Name	Anti-ErbB (HER3) Rabbit Monoclonal Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-ErbB (HER3) Rabbit Monoclonal Antibody catalog # M00539. Tested in WB, IP applications. This antibody reacts with Human, Rat.
Application	IP, WB
Clonality	Monoclonal BIB-5
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P21860

### Technical Details

Immunogen	A synthesized peptide derived from human ErbB (HER3)
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IP 1:30

## Anti-ErbB (HER3) Rabbit Monoclonal Antibody (M00539) Images



Western blot analysis of ErbB3 (HER3) expression in MCF-7 cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ErbB (HER3) Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.