

Anti-Cyclin E1 CCNE1 Rabbit Monoclonal Antibody

Catalog Number: M00543-1


Overview

Product Name	Anti-Cyclin E1 CCNE1 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Cyclin E1 CCNE1 Rabbit Monoclonal Antibody catalog # M00543-1. Tested in WB, ICC/IF, IP, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IP, IF, ICC, WB
Clonality	Monoclonal DDF-3
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P24864

Technical Details

Immunogen	A synthesized peptide derived from human Cyclin E1
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 ICC/IF 1:50-1:200 IP 1:50 FC 1:50

Anti-Cyclin E1 CCNE1 Rabbit Monoclonal Antibody (M00543-1) Images



Western blot analysis of Cyclin E1 expression in HeLa cell lysates.

1 Publications Citing This Product

1. PubMed ID: 23145020, Wu Y, Tao Y, Chen Y, Xu W. Plos One. 2012;7(11):E48917. Doi: 10.1371/Journal.Pone.0048917. Epub 2012 Nov 7. Rhoc Regulates The Proliferation Of Gastric Cancer Cells Through Interaction With Iqgap1.

Visit bosterbio.com/anti-cyclin-e1-rabbit-monoclonal-antibody-m00543-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cyclin E1 CCNE1 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.