

Anti-Wnt5a Rabbit Monoclonal Antibody

Catalog Number: M00549


Overview

Product Name	Anti-Wnt5a Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Wnt5a Rabbit Monoclonal Antibody catalog # M00549. Tested in WB, ICC/IF, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, ICC, WB
Clonality	Monoclonal ABEE-23
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P41221


Technical Details

Immunogen	A synthesized peptide derived from human Wnt5a
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000 ICC/IF 1:50-200 FC 1:20


Anti-Wnt5a Rabbit Monoclonal Antibody (M00549) Images




All lanes use the Antibody at 1:3K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:3K dilution for 1 hour at room temperature.





The mechanism of regulatory relationships among EZH2, SFRP1, and WNT5A. (A) The miR26a level in C4-2B and C4-2B EnzR cells. n = 3, mean \pm SD, *** p



All lanes use the Antibody at 1:3K dilution for 1 hour at room temperature.

Western blot analysis of Wnt5a expression in (1) HeLa cell lysate; (2) SK-OV-3 cell lysate.



All lanes use the Antibody at 1:3K dilution for 1 hour at room temperature.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Wnt5a Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.