

## Anti-PTPRC / CD45 Reference Antibody (Novartis patent anti-CD45)

Catalog Number: M00555-8


### Overview

Product Name	Anti-PTPRC / CD45 Reference Antibody (Novartis patent anti-CD45)
Reactive Species	Human
Description	Boster Bio Anti-PTPRC / CD45 Reference Antibody (Novartis patent anti-CD45) (Catalog # M00555-8). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.412EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P08575


### Technical Details

Isotype	IgG1
Form	Liquid


## Anti-PTPRC / CD45 Reference Antibody (Novartis patent anti-CD45) (M00555-8) Images



Anti-PTPRC/CD45 Reference Antibody (Novartis patent anti-CD45) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-PTPRC/CD45 Reference Antibody (Novartis patent anti-CD45) is more than 95%



Immobilized human IL-18 His at 4 µg/mL can bind Anti-PTPRC/CD45 Reference Antibody (Novartis patent anti-CD45)

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PTPRC / CD45 Reference Antibody (Novartis patent anti-CD45)

For Research Use Only. Not for use in diagnostic procedures.