

## Anti-EPHA2 Antibody (Monoclonal, 33E23)

Catalog Number: M00578-3

### About EPHA2

EPHA2(efhrin type-A receptor 2) also known as ECK, is a protein that in humans is encoded by the EPHA2 gene. This gene belongs to the ephrin receptor subfamily of the protein-tyrosine kinase family. Receptors in the EPH subfamily typically have a single kinase domain and an extracellular region containing a Cys-rich domain and 2 fibronectin type III repeats. By somatic cell hybrid analysis and fluorescence in situ hybridization, the EPHA2 gene is mapped to chromosome 1p36.1. EPHA2 was readily detectable in human lens fiber cells using immunoblot and immunohistochemistry. EGFR and EPHA2 mediated HCV entry by regulating CD81 -claudin-1 (CLDN1) coreceptor associations and viral glycoprotein-dependent membrane fusion.

### Overview

Product Name	Anti-EPHA2 Antibody (Monoclonal, 33E23)
Reactive Species	Human
Description	Boster Bio Anti-EPHA2 Antibody (Monoclonal, 33E23) catalog # M00578-3. Tested in WB, ICC/IF, IP, Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IP, IF, ICC, WB
Clonality	Monoclonal 33E23
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P29317

### Technical Details

Immunogen	Recombinant protein within human EPHA2 aa 12-548.
Form	Liquid
Concentration	500 ug/ml
Purification	Protein A affinity purified.

**Suggested Dilutions**

Western blot, 1:500-2000  
Immunocytochemistry/Immunofluorescence, 1:50-200  
ImmunoPrecipitation, 1:50  
Flow Cytometry (Fixed), 1:50-200

**Submit a product review to Biocompare.com**

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EPHA2 Antibody (Monoclonal, 33E23)

For Research Use Only. Not for use in diagnostic procedures.