

Anti-CCL5 / RANTES Reference Antibody (NI-0701)

Catalog Number: M00617-3


Overview

Product Name	Anti-CCL5 / RANTES Reference Antibody (NI-0701)
Reactive Species	Cynomolgus monkey, Human, Mouse
Description	Boster Bio Anti-CCL5 / RANTES Reference Antibody (NI-0701) (Catalog # M00617-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human, Mouse. Endotoxin: < 0.424EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in 100mM Pro-Ac, 20mM Arg, pH5.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P13501


Technical Details

Isotype	IgG2SA
Form	Liquid

Anti-CCL5 / RANTES Reference Antibody (NI-0701) (M00617-3) Images



Anti-CCL5/RANTES Reference Antibody (NI-0701) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-CCL5/RANTES Reference Antibody (NI-0701) is more than 98.94%

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CCL5 / RANTES Reference Antibody (NI-0701)

For Research Use Only. Not for use in diagnostic procedures.