

Anti-M-CSF (CSF1) Mouse Monoclonal Antibody [Clone ID: OTI2D6]

Catalog Number: M00620-1


Overview

Product Name	Anti-M-CSF (CSF1) Mouse Monoclonal Antibody [Clone ID: OTI2D6]
Reactive Species	Human
Description	Boster Bio CSF1 mouse monoclonal antibody, clone OTI2D6 (formerly 2D6). Catalog# M00620-1. Tested in IHC, WB. This antibody reacts with Human.
Conjugate	Unconjugated
Application	IHC, WB
Clonality	Monoclonal OTI2D6
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P09603


Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 223-496 of human CSF1 (NP_000748) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IHC 1:150


Anti-M-CSF (CSF1) Mouse Monoclonal Antibody [Clone ID: OTI2D6] (M00620-1) Images



Immunohistochemical staining of paraffin-embedded Carcinoma of Human lung tissue using anti-CSF1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)



Immunohistochemical staining of paraffin-embedded Human tonsil within the normal limits using anti-CSF1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris)



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CSF1 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CSF1.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-M-CSF (CSF1) Mouse Monoclonal Antibody [Clone ID: OTI2D6]

For Research Use Only. Not for use in diagnostic procedures.