

Anti-CSF1 / M-CSF Reference Antibody (lacnotuzumab)

Catalog Number: M00620-2


Overview

Product Name	Anti-CSF1 / M-CSF Reference Antibody (lacnotuzumab)
Reactive Species	Human
Description	Boster Bio Anti-CSF1 / M-CSF Reference Antibody (lacnotuzumab) (Catalog # M00620-2). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.426EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P09603


Technical Details

Isotype	IgG1
Form	Liquid

Anti-CSF1 / M-CSF Reference Antibody (Iacnotuzumab) (M00620-2) Images



Anti-CSF1/M-CSF Reference Antibody (Iacnotuzumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-CSF1/M-CSF Reference Antibody (Iacnotuzumab) is more than 95%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CSF1 / M-CSF Reference Antibody (Iacnotuzumab)

For Research Use Only. Not for use in diagnostic procedures.