

Anti-Factor X (F10) Mouse Monoclonal Antibody [Clone ID: OTI1D12]

Catalog Number: M00641


Overview

Product Name	Anti-Factor X (F10) Mouse Monoclonal Antibody [Clone ID: OTI1D12]
Reactive Species	Human
Description	Boster Bio F10 mouse monoclonal antibody, clone OTI1D12. Catalog# M00641. Tested in WB. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal OTI1D12
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C as received.
Host	Mouse
Uniprot ID	P00742

Technical Details

Immunogen	Human recombinant protein fragment corresponding to amino acids 41-488 of human F10 (NP_000495) produced in E.coli.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB: 1:2000

Anti-Factor X (F10) Mouse Monoclonal Antibody [Clone ID: OTI1D12] (M00641) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY F10 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-F10 (1:2000).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Factor X (F10) Mouse Monoclonal Antibody [Clone ID: OTI1D12]

For Research Use Only. Not for use in diagnostic procedures.