

## Anti-HMGCR/Hmg Coa Reductase Rabbit Monoclonal Antibody

Catalog Number: M00643


### Overview

Product Name	Anti-HMGCR/Hmg Coa Reductase Rabbit Monoclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HMGCR/Hmg Coa Reductase Rabbit Monoclonal Antibody catalog # M00643. Tested in WB, IP applications. This antibody reacts with Human, Mouse, Rat.
Application	IP, WB
Clonality	Monoclonal AOCC-8
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P04035

### Technical Details

Immunogen	A synthesized peptide derived from human HMGCR
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-1:2000 IP 1:50

## Anti-HMGCR/Hmg Coa Reductase Rabbit Monoclonal Antibody (M00643) Images



Western blot analysis of HMGCR expression in Jurkat cell lysate.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HMGCR/Hmg Coa Reductase Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.