

## Anti-TIM-3 / HAVCR2 / CD366 Reference Antibody (sabatolimab)

Catalog Number: M00657-3


### Overview

Product Name	Anti-TIM-3 / HAVCR2 / CD366 Reference Antibody (sabatolimab)
Reactive Species	Cynomolgus monkey, Human
Description	Boster Bio Anti-TIM-3 / HAVCR2 / CD366 Reference Antibody (sabatolimab) (Catalog # M00657-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.434EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	Q8TDQ0


### Technical Details

Isotype	IgG4SP
Form	Liquid

## Anti-TIM-3 / HAVCR2 / CD366 Reference Antibody (sabatolimab) (M00657-3) Images



Anti-TIM-3/HAVCR2/CD366 Reference Antibody (sabatolimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TIM-3/HAVCR2/CD366 Reference Antibody (sabatolimab) is more than 98.49%

**Submit a product review to Biocompare.com**

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TIM-3 / HAVCR2 / CD366 Reference Antibody (sabatolimab)

For Research Use Only. Not for use in diagnostic procedures.