

# Anti-Human Emerin DyLight® 488 conjugated EMD Antibody(monoclonal, 5A10)

Catalog Number: M00714-Dyl488

## Overview

Product Name	Anti-Human Emerin DyLight® 488 conjugated EMD Antibody(monoclonal, 5A10)
Reactive Species	Human
Description	Boster Bio Anti-Human Emerin DyLight® 488 conjugated EMD Antibody (monoclonal, 5A10) catalog # M00714-Dyl488. Tested in Flow Cytometry applications. This antibody reacts with Human.
Conjugate	DyLight®488
Application	Flow Cytometry
Clonality	Monoclonal 5A10
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	P50402

## Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human Emerin, different from the related mouse sequence by eight amino acids, and from the related rat sequence by nine amino acids.
Predicted Reactive Species	Human
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG1
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used:

Flow Cytometry, 1-3ug/1x10<sup>6</sup> cells

## Anti-Human Emerin DyLight® 488 conjugated EMD Antibody(monoclonal, 5A10) (M00714-Dyl488) Images

---



Boster Kit Box

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Human Emerin DyLight 488 conjugated EMD Antibody(monoclonal, 5A10)