

Anti-SNAI 1 Monoclonal Antibody

Catalog Number: M00716


Overview

Product Name	Anti-SNAI 1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-SNAI 1 Monoclonal Antibody catalog # M00716. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 6D2
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O95863


Technical Details

Immunogen	Purified recombinant fragment of human SNAI 1 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Affinity purification
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-SNAI 1 Monoclonal Antibody (M00716) Images



Western blotting validation for Anti-SNAI 1 Monoclonal Antibody M00716 Western Blot (WB) analysis using SNAI 1 Monoclonal Antibody against NTERA-2 cell lysate. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



ELISA validation of SNAI1 using Anti-SNAI 1 Monoclonal Antibody (M00716). ELISA analysis of SNAI 1 antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SNAI 1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.