

## Anti-Apolipoprotein A I APOA1 Antibody Picoband® (monoclonal, 17G5) Fluoro647 Conjugated

Catalog Number: M00717-1-Fluoro647

### About APOA1

Apolipoprotein A-1, also known as APOA1, is a human protein with a specific role in lipid metabolism. It binds to lipopolysaccharide or endotoxin, and has a major role in the anti-endotoxin function of HDL. The gene is mapped to 11q23. And it is a single polypeptide chain with 243 amino acid residues of known primary amino acid sequence. The ApoA-I protein promotes cholesterol efflux from tissues to the liver for excretion. It is a cofactor for lecithin cholesterolacyltransferase (LCAT) which is responsible for the formation of most plasma cholesteryl esters. ApoA-I is also isolated as a prostacyclin (PGI<sub>2</sub>) stabilizing factor, and thus may have an anticlotting effect. Defects in the gene encoding it are associated with HDL deficiencies, including Tangier disease, and with systemic non-neuropathic amyloidosis. Additionally, ApoA-I overexpression promotes macrophage-specific reverse cholesterol transport.

### Overview

Product Name	Anti-Apolipoprotein A I APOA1 Antibody Picoband® (monoclonal, 17G5) Fluoro647 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (WB, ELISA (Cap)). Customers may select suitable applications according to their experimental needs.
Clonality	Monoclonal 17G5
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Mouse
Uniprot ID	P02647

### Technical Details

Immunogen	E. coli-derived mouse Apolipoprotein A I recombinant protein (Position: D25-Q264). Mouse Apolipoprotein A I shares 64% and 68.6% amino acid (aa) sequence identity with human and rat Apolipoprotein A I, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG1
Form	Liquid
Concentration	0.5 mg/mL
Purification	Affinity-chromatography

Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 34103620, Xiang Y,Zheng Y,Liu S,Liu G,Li Z,Dong W.Comparison of the sensitivity of Western blotting between PVDF and NC membranes.Sci Rep.2021 Jun 8;11(1):12022.doi:10.1038/s41598-021-91521-8.PMID:34103620;PMCID:PMC8187435.

Visit [bosterbio.com/anti-apolipoprotein-picoband-trade-antibody-monoclonal-m00717-1-boster.html](https://bosterbio.com/anti-apolipoprotein-picoband-trade-antibody-monoclonal-m00717-1-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Apolipoprotein A I APOA1 Antibody ® (monoclonal, 17G5) - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.