

Anti-STAT4 Mouse Monoclonal Antibody [Clone ID: OTI2F3]

Catalog Number: M00734-1


Overview

Product Name	Anti-STAT4 Mouse Monoclonal Antibody [Clone ID: OTI2F3]
Reactive Species	Human, Mouse, Rat
Description	Boster Bio STAT4 mouse monoclonal antibody, clone OTI2F3 (formerly 2F3). Catalog# M00734-1. Tested in FC, IF, IHC, WB. This antibody reacts with Human, Mouse, Rat.
Conjugate	Unconjugated
Application	Flow Cytometry, IF, IHC, WB
Clonality	Monoclonal OTI2F3
Formulation	PBS (pH 7.3) containing 1% stabilizing protein, 50% glycerol and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q14765


Technical Details

Immunogen	Full length human recombinant protein of human STAT4 (NP_003142) produced in HEK293T cell.
Isotype	IgG1
Concentration	1 mg/ml
Purification	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Suggested Dilutions	WB 1:2000 IHC 1:150 IF 1:100 Flow Cytometry 1:100


Anti-STAT4 Mouse Monoclonal Antibody [Clone ID: OTI2F3] (M00734-1) Images




HEK293T cells transfected with either STAT4 (Myc-DDK-tagged) overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-STAT4 antibody (M00734-1)




Anti-STAT4 mouse monoclonal antibody (M00734-1) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY STAT4.



Immunohistochemical staining of paraffin-embedded Human lymphoma tissue using anti-STAT4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



Immunohistochemical staining of paraffin-embedded Human tonsil within the normal limits using anti-STAT4 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer)



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY STAT4 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-STAT4.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-STAT4 Mouse Monoclonal Antibody [Clone ID: OT12F3]

For Research Use Only. Not for use in diagnostic procedures.