

Anti-hCG (Human Chorionic Gonadotropin) CGA Rabbit Monoclonal Antibody, Clone#RM330

Catalog Number: M00739-2


Overview

Product Name	Anti-hCG (Human Chorionic Gonadotropin) CGA Rabbit Monoclonal Antibody, Clone#RM330
Reactive Species	Human
Description	Boster Bio Anti-hCG (Human Chorionic Gonadotropin) CGA Rabbit Monoclonal Antibody, Clone#RM330 (Catalog # M00739-2). Tested in WB, IHC, ELISA applications. This antibody reacts with Human.
Application	ELISA, IHC, WB
Clonality	Monoclonal RM330
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01215


Technical Details

Immunogen	Native hCG protein purified from Human Urine.
Cross Reactivity	This antibody reacts to hCG (human Chorionic Gonadotropin). The antibody recognizes both the native complex and the beta subunit of hCG.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	Immunohistochemistry (IHC): 1:500-1:1000 dilution WB: 1:5000-1:20000 dilution ELISA: 1:1000-1:10000.

Anti-hCG (Human Chorionic Gonadotropin) CGA Rabbit Monoclonal Antibody, Clone#RM330 (M00739-2) Images



Western Blotting result Western Blot of hCG protein purified from human urine using anti-hCG rabbit monoclonal antibody (Clone RM330) at a 1:20,000 dilution.



IHC result Immunohistochemical staining of formalin fixed and paraffin embedded human placenta tissue section using anti-hCG rabbit monoclonal antibody (Clone RM330) at a 1:1000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-hCG (Human Chorionic Gonadotropin) CGA Rabbit Monoclonal Antibody, Clone#RM330

For Research Use Only. Not for use in diagnostic procedures.