

Anti-TNFSF1 / TNFb/ LT alpha Reference Antibody (pateclizumab)

Catalog Number: M00746-3


Overview

| | |
|----------------------|---|
| Product Name | Anti-TNFSF1 / TNFb/ LT alpha Reference Antibody (pateclizumab) |
| Reactive Species | Cynomolgus monkey, Human |
| Description | Boster Bio Anti-TNFSF1 / TNFb/ LT alpha Reference Antibody (pateclizumab) (Catalog # M00746-3). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Cynomolgus monkey, Human. Endotoxin: < 0.444EU/ug,determined by LAL method. Expression system: CHO Cell |
| Application | ELISA, Flow Cytometry, Functional Assay, Kinetics |
| Clonality | Monoclonal |
| Formulation | Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column. |
| Storage Instructions | Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles. |
| Uniprot ID | P01374 |


Technical Details

| | |
|---------|--------|
| Isotype | IgG1 |
| Form | Liquid |

Anti-TNFSF1 / TNFb/ LT alpha Reference Antibody (pateclizumab) (M00746-3) Images



Anti-TNFSF1/TNFb/ LT alpha Reference Antibody (pateclizumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-TNFSF1/TNFb/ LT alpha Reference Antibody (pateclizumab) is more than 95%

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TNFSF1 / TNFb/ LT alpha Reference Antibody (pateclizumab)

For Research Use Only. Not for use in diagnostic procedures.