

Anti-CCL20 Reference Antibody (GSK3050002)

Catalog Number: M00748

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-CCL20 Reference Antibody (GSK3050002)
Reactive Species	Human
Description	Boster Bio Anti-CCL20 Reference Antibody (GSK3050002) (Catalog # M00748). Tested in Flow Cytometry, ELISA, FTA. This antibody reacts with Human. Endotoxin: < 0.445EU/ug,determined by LAL method. Expression system: CHO Cell
Application	ELISA, Flow Cytometry, Functional Assay, Kinetics
Clonality	Monoclonal
Formulation	Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.
Storage Instructions	Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.
Uniprot ID	P78556


Technical Details

Isotype	IgG4
Form	Liquid

Anti-CCL20 Reference Antibody (GSK3050002) (M00748) Images



Anti-CCL20 Reference Antibody (GSK3050002) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-CCL20 Reference Antibody (GSK3050002) is more than 98.28%

[Anti-CCL20 Reference Antibody \(GSK3050002\)](#)

For Research Use Only. Not for use in diagnostic procedures.