

Anti-ZAP70 Antibody Picoband® (monoclonal, 9D5) Biotin Conjugated

Catalog Number: M00754-5-Biotin

About ZAP70

ZAP-70 (Zeta-chain-associated protein kinase 70) encodes an enzyme belonging to the protein tyrosine kinase family, and it plays a role in T-cell development and lymphocyte activation. This enzyme, which is phosphorylated on tyrosine residues upon T-cell antigen receptor (TCR) stimulation, functions in the initial step of TCR-mediated signal transduction in combination with the Src family kinases, Lck and Fyn. This enzyme is also essential for thymocyte development. Mutations in this gene cause selective T-cell defect, a severe combined immunodeficiency disease characterized by a selective absence of CD8-positive T-cells. Two transcript variants that encode different isoforms have been found for this gene.

Overview

Product Name	Anti-ZAP70 Antibody Picoband® (monoclonal, 9D5) Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Monoclonal 9D5
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Mouse
Uniprot ID	P43403

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human ZAP70, different from the related mouse sequence by two amino acids.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Mouse IgG2b
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZAP70 Antibody (monoclonal, 9D5) - Biotin

For Research Use Only. Not for use in diagnostic procedures.