

## Anti-TGF beta Receptor II TGFBR2 Rabbit Monoclonal Antibody

Catalog Number: M00759


### Overview

Product Name	Anti-TGF beta Receptor II TGFBR2 Rabbit Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-TGF beta Receptor II TGFBR2 Rabbit Monoclonal Antibody catalog # M00759. Tested in WB application. This antibody reacts with Human.
Application	WB
Clonality	Monoclonal ABFE-20
Formulation	Rabbit IgG in stabilizing components, phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P37173


### Technical Details

Immunogen	A synthesized peptide derived from human TGF beta Receptor II
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Suggested Dilutions	WB 1:500-2000


## Anti-TGF beta Receptor II TGFBR2 Rabbit Monoclonal Antibody (M00759) Images



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of TGF beta Receptor II expression in A549 cell lysate.

### 1 Publications Citing This Product

1. PubMed ID: 27019660, Effect of Kujjie Granule on the Expression of TGF- $\beta$ /Smads Signaling Pathway in Patients with Ulcerative Colitis

Visit [bosterbio.com/anti-tgf-beta-receptor-ii-rabbit-monoclonal-antibody-m00759-boster.html](http://bosterbio.com/anti-tgf-beta-receptor-ii-rabbit-monoclonal-antibody-m00759-boster.html) to see all 1 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TGF beta Receptor II TGFBR2 Rabbit Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.