

## Anti-Aurora-B (AURKB) Rabbit Monoclonal Antibody, Clone#RM278

Catalog Number: M00762-4


### Overview

Product Name	Anti-Aurora-B (AURKB) Rabbit Monoclonal Antibody, Clone#RM278
Reactive Species	Human
Description	Boster Bio Anti-Aurora-B (AURKB) Rabbit Monoclonal Antibody, Clone#RM278 (Catalog # M00762-4). Tested in IHC, WB applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Monoclonal RM278
Formulation	50% Glycerol/PBS with 1% stabilizing protein and 0.09% sodium azide This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96GD4


### Technical Details

Immunogen	A peptide corresponding to the N-terminus of human Aurora kinase B.
Cross Reactivity	This antibody reacts to human Aurora kinase B.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	Protein A affinity purified from an animal origin-free culture supernatant
Suggested Dilutions	Immunohistochemistry (IHC): 1:200-1:500 dilution WB: 1:1000-1:2000 dilution.

## Anti-Aurora-B (AURKB) Rabbit Monoclonal Antibody, Clone#RM278 (M00762-4) Images



Western Blotting result Western Blot of HeLa cells nontreated (-) or treated (+) with Nocodazole (NOC), using Anti-Aurora B Rabbit Monoclonal Antibody (Clone RM278) at a 1:1000 dilution.



IHC result Immunohistochemical staining of formalin fixed and paraffin embedded human tonsil tissue sections using Anti-Aurora B Rabbit Monoclonal Antibody (Clone RM278) at a 1:500 dilution.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Aurora-B (AURKB) Rabbit Monoclonal Antibody, Clone#RM278

For Research Use Only. Not for use in diagnostic procedures.